Caretray Trolleys & Racks

s 34	Stainless Steel, Low Level Trolleys
s 36	Stainless Steel, Standard Level Trolleys
s 38	Stainless Steel, High Level Trolleys
s 40	Stainless Steel, Mobile, Racks
s 42	Stainless Steel, Static, Racks
s 44	Mild Steel, Low Level Trolleys
s 46	Mild Steel, Standard Level Trolleys
s 48	Mild Steel, High Level Trolleys
s 50	Mild Steel, Mobile Racks
s 52	Mild Steel, Static Racks
s 54	Accessories

Caretray Trolleys - Stainless Steel, Low Level, Single & Double Column

Features

- Low level, single & double column easy clean procedure trolleys
- · Quick & easy disassembly, assists thorough & deep cleaning
- Seamlessly welded, one-piece 304 grade stainless steel frame
- Recyclable polypropylene one-piece top & side panels, smooth surfaces & radiused edges assists cleaning (compatible with steam cleaning)
 - Top Panel Raised edges on three sides provide defined grab points to act as push handles, whilst an integral rim at the front contains spillages (max load 30kg)
 - Side Panels Incorporate runners with positive stops preventing trays from being accidentally withdrawn (symmetrical design allows tray insertion from either side)
- Carerail fitted to either side support a wide range of accessories
- Translucent, impact-resistant trays & dividers, smooth surfaces & radiused edges assists cleaning (compatible with steam cleaning / high temperature machine washing)
- Trays feature an integral handle with a plastic clip-on label cover (coloured labels supplied - white, blue, green, yellow & red)
- Each tray supplied with 2 x long & 2 x short dividers
- · Max load per tray 5kg
- Polyurethane ring buffers offer protection to both the trolley & fabric of the building, smooth surfaces assist cleaning
- · 100mm pressed steel castors, 2 off braking

Dimensions (w x d x h)

- Overall Dimensions
 - Single Column 550 x 550 x 807mm
 - Double Column 1005 x 550 x 807mm
- Plastic Top
 - Single Column 455 x 455 (internal) at height 782mm
 - Double Column 915 x 455 (internal) at height 782mm
- Travs
 - Shallow Tray 355 x 405 x 100mm (internal)
 - Deep Tray 355 x 405 x 150mm (internal)

Item	Description	Weight
CTS1/08/1S2D	1 x 100mm, 2 x 150mm Trays	21.2kg
CTS1/08/3D	3 x 150mm Trays	21.5kg
CTS1/08/3S1D	3 x 100mm, 1 x 150mm Trays	22.3kg
CTS1/08/4S	4 x 100mm Trays	21.9kg



CTS1/08/3S1D - shown with optional CT/WF, CT/ALC & GMBH/15

Easy disassembly & cleaning

Options / Accessories

- Locking & Security page 54
- Overbridge & Tilt Bins page 60
- Defibrillator Platform page 61
- Carerail Mounted page 62
- Laptop Holders page 68
- Castors & Ring Buffers page 70
 Trays, Shelves, Dividers & Labels page 71







Caretray Trolleys - Stainless Steel, Standard Level, Single & Double Column

Features

- Standard level, single & double column easy clean procedure trollevs
- · Quick & easy disassembly, assists thorough & deep cleaning
- · Seamlessly welded, one-piece 304 grade stainless steel frame
- Recyclable polypropylene one-piece top & side panels, smooth surfaces & radiused edges assists cleaning (compatible with steam cleaning)
 - Top Panel Raised edges on three sides provide defined grab points to act as push handles, whilst an integral rim at the front contains spillages (max load 30kg)
 - Side Panels Incorporate runners with positive stops preventing trays from being accidentally withdrawn (symmetrical design allows tray insertion from either side)
- Carerails fitted to either side support a wide range of accessories
- Translucent, impact-resistant trays & dividers, smooth surfaces radiused edges assists cleaning (compatible with steam cleaning / high temperature machine washing)
- Trays feature an integral handle with a plastic clip-on label cover (coloured labels supplied white, blue, green, yellow & red)
- Each tray supplied with 2 x long & 2 x short dividers
- Max load per tray 5kg
- Polyurethane ring buffers offer protection to both the trolley & fabric of the building
- · 100mm pressed steel castors, 2 off braking

Dimensions (w x d x h)

- Overall Dimensions
 - Single Column 550 x 550 x 993mm
 - Double Column 1005 x 550 x 993mm
- Plastic Top
 - Single Column 455 x 455 (internal) at height 968mm
 - Double Column 915 x 455 (internal) at height 968mm
- - Shallow Tray 355 x 405 x 100mm (internal)
 - Deep Tray 355 x 405 x 150mm (internal)










CTS1/10/3S2D



Options / Accessories • Locking & Security - page 54 • Overbridge & Tilt Bins - page 60 • Defibrillator Platform - page 61

- Carerail Mounted page 62

- Laptop Holders page 68
 Castors & Ring Buffers page 70
 Trays, Shelves, Dividers & Labels page 71



Item	Description	Weight
CTS2/10/3S6D	3 x 100mm, 6 x 150mm Trays	46.6kg
CTS2/10/4S5D	4 x 100mm, 5 x 150mm Trays	46.2kg
CTS2/10/5S4D	5 x 100mm, 4 x 150mm Trays	45.9kg
CTS2/10/6S4D	6 x 100mm, 4 x 150mm Trays	47.3kg
CTS2/10/8D	8 x 150mm Trays	45.8kg
CTS2/10/8S2D	8 x 100mm, 2 x 150mm Trays	46.6kg



CT/IV, CT/CR2, GMBH/15 & CT/ALC

Caretray Trolleys - Stainless Steel, High Level, Single & Double Column

Features

- High level, single & double column easy clean procedure trolleys
- · Quick & easy disassembly, assists thorough & deep cleaning
- Seamlessly welded, one-piece 304 grade stainless steel frame
- Recyclable polypropylene one-piece top & side panels, smooth surfaces & radiused edges assists cleaning (compatible with steam cleaning)
 - Top Panel Raised edges on three sides provide defined grab points to act as push handles, whilst an integral rim at the front contains spillages (max load 30kg)
 - Side Panels Incorporate runners with positive stops preventing trays from being accidentally withdrawn (symmetrical design allows tray insertion from either side)
- Carerail fitted to either side support a wide range of accessories
- Translucent, impact-resistant trays & dividers, smooth surfaces & radiused edges assists cleaning (compatible with steam cleaning / high temperature machine washing)
- Trays feature an integral handle with a plastic clip-on label cover (coloured labels supplied white, blue, green, yellow & red)
- Each tray supplied with 2 x long & 2 x short dividers
- Max load per tray 5kg
- Polyurethane ring buffers offer protection to both the trolley & fabric of the building
- 100mm pressed steel castors, 2 off braking

Dimensions (w x d x h)

- Overall Dimensions
 - Single Column 550 x 550 x 1410mm
 - Double Column 1005 x 550 x 1410mm
- Plastic Top
 - Single Column 455 x 455 (internal) at height 1385mm
 - Double Column 915 x 455 (internal) at height 1385mm
- Trays
 - Shallow Tray 355 x 405 x 100mm (internal)
 - Deep Tray 355 x 405 x 150mm (internal)



CTS1/14/4S3D - shown with optional CT/DH/A4

Easy disassembly & cleaning



CTS1/14/4S3D

Item	Description	Weight
CTS1/14/2S5D	2 x 100mm, 5 x 150mm Trays	36.2kg
CTS1/14/3S4D	3 x 100mm, 4 x 150mm Trays	35.8kg
CTS1/14/4S3D	4 x 100mm, 3 x 150mm Trays	35.4kg
CTS1/14/5S3D	5 x 100mm, 3 x 150mm Trays	37.8kg
CTS1/14/6D	6 x 150mm Trays	35.0kg
CTS1/14/8S1D	8 x 100mm 1 x 150mm Travs	37 8kg

Options / Accessories

- Locking & Security page 54Carerail Mounted page 62

- Castors & Ring Buffers page 70
 Trays, Shelves, Dividers & Labels page 71



Item Description CTS2/14/4S10D 2 x 100mm, 10 x 150mm Trays CTS2/14/6S8D 6 x 100mm, 8 x 150mm Trays CTS2/14/8S6D 8 x 100mm, 6 x 150mm Trays CTS2/14/10S6D 10 x 100mm, 6 x 150mm Trays CTS2/14/12D 12 x 150mm Trays CTS2/14/16S2D 16 x 100mm, 2 x 150mm Trays 66.6kg

Caretray Racks - Stainless Steel, Mobile, Single & Double Column

Features

- Mobile, single & double column easy clean racks
- · Quick & easy disassembly, assists thorough & deep cleaning
- Seamlessly welded, one-piece 304 grade stainless steel frame
- Recyclable polypropylene one-piece top & side panels, smooth surfaces & radiused edges assists cleaning (compatible with steam cleaning)
 - Top Panel Raised edges on three sides provide defined grab points to act as push handles, whilst an integral rim at the front contains spillages (max load 30kg)
 - Side Panels Incorporate runners with positive stops preventing trays from being accidentally withdrawn (symmetrical design allows tray insertion from either side)
- Translucent, impact-resistant travs & dividers, smooth surfaces & radiused edges assists cleaning (compatible with steam cleaning / high temperature machine washing)
- Trays feature an integral handle with a plastic clip-on label cover (coloured labels supplied white, blue, green, yellow & red)
- Each tray supplied with 2 x long & 2 x short dividers
- Max load per trav 5kg
- Polyurethane ring buffers offer protection to both the rack & fabric of the building
- 100mm pressed steel castors. 2 off braking

Dimensions (w x d x h)

- Overall Dimensions
 - Single Column 550 x 550 x 1782mm
 - Double Column 1005 x 550 x 1782mm
- Plastic Top
 - Single Column 455 x 455 (internal) at height 1757mm
 - Double Column 915 x 455 (internal) at height 1757mm
- Trays

Item

CRSM1/18/2S6D

CRSM1/18/3S6D

CRSM1/18/4S5D

CRSM1/18/5S5D

CRSM1/18/6S4D

CRSM1/18/8D

• Shallow Tray - 355 x 405 x 100mm (internal)

Description

8 x 150mm Trays

2 x 100mm, 6 x 150mm Trays

3 x 100mm, 6 x 150mm Trays

4 x 100mm, 5 x 150mm Trays

5 x 100mm, 5 x 150mm Trays

6 x 100mm, 4 x 150mm Trays

• Deep Tray - 355 x 405 x 150mm (internal)



Weight

41.7kg

44.1kg

42.8kg

44.3kg

43.9kg

43.3kg







CRSM1/18/6S4D

Options / Accessories

- Castors & Ring Buffers page 70Trays, Shelves, Dividers & Labels page 71



Item Description Weight CRSM2/18/4S12D 4 x 100mm, 12 x 150mm Trays 72.1kg CRSM2/18/6S12D 6 x 100mm, 12 x 150mm Trays 77.1kg CRSM2/18/8S10D 8 x 100mm, 10 x 150mm Trays 75.4kg CRSM2/18/10S10D 10 x 100mm, 10 x 150mm Trays 78.4kg CRSM2/18/12S8D 12 x 100mm, 8 x 150mm Trays 77.7kg CRSM2/18/16D 16 x 150mm Trays 74.6kg



Caretray Racks - Stainless Steel, Static, Single & Double Column

Weiaht

40.0kg

Features

- Static, single & double column easy clean racks
- · Quick & easy disassembly, assists thorough & deep cleaning
- · Seamlessly welded, one-piece 304 grade stainless steel frame
- Recyclable polypropylene one-piece top & side panels, smooth surfaces & radiused edges assists cleaning (compatible with steam cleaning)
 - Side Panels Incorporate runners with positive stops preventing trays from being accidentally withdrawn (symmetrical design allows tray insertion from either side)
- Translucent, impact-resistant trays & dividers, smooth surfaces & radiused edges assists cleaning (compatible with steam cleaning / high temperature machine washing)
- Trays feature an integral handle with a plastic clip-on label cover (coloured labels supplied white, blue, green, yellow & red)
- Each tray supplied with 2 x long & 2 x short dividers
- · Max load per tray 5kg
- Adjustable feet allow levelling on uneven floors

Dimensions (w x d x h)

Overall Dimensions

CDCC4/4C/0CCD

CRSS1/16/8D

- Single Column 550 x 550 x 1685mm
- Double Column 1005 x 550 x 1685mm
- Plastic Top
 - Single Column 455 x 455 (internal) at height 1645mm
 - Double Column 915 x 455 (internal) at height 1645mm
- Travs
 - Shallow Tray 355 x 405 x 100mm (internal)
 - Deep Tray 355 x 405 x 150mm (internal)



CRSS 1/10/250D	2 x 100mm, 6 x 150mm mays	36.4Kg
CRSS1/16/3S6D	3 x 100mm, 6 x 150mm Trays	40.8kg
CRSS1/16/4S5D	4 x 100mm, 5 x 150mm Trays	39.5kg
CRSS1/16/5S5D	5 x 100mm, 5 x 150mm Trays	41.0kg
CRSS1/16/6S4D	6 x 100mm, 4 x 150mm Trays	40.6kg

8 x 150mm Trays

0 . . 400---- C . . 450--- T-- . .

Description

Easy disassembly & cleaning

Options / Accessories
Trays, Shelves, Dividers & Labels - page 71



Item	Description	Weight
CRSS2/16/4S12D	4 x 100mm, 12 x 150mm Trays	71.1kg
CRSS2/16/6S12D	6 x 100mm, 12 x 150mm Trays	74.0kg
CRSS2/16/8S10D	8 x 100mm, 10 x 150mm Trays	72.3kg
CRSS2/16/10S10D	10 x 100mm, 10 x 150mm Trays	75.3kg
CRSS2/16/12S8D	12 x 100mm, 8 x 150mm Trays	74.6kg
CRSS2/16/16D	16 x 150mm Trays	71.5kg



Caretray Trolleys - Mild Steel, Low Level, Single & Double Column

Features

- Low level, single & double column easy clean procedure trolleys
- · Quick & easy disassembly, assists thorough & deep cleaning
- Recyclable polypropylene one-piece top & side panels, smooth surfaces & radiused edges assists cleaning (compatible with steam cleaning)
 - Top Panel Raised edges on three sides provide defined grab points to act as push handles, whilst an integral rim at the front contains spillages (max load 30kg)
 - Side Panels Incorporate runners with positive stops preventing trays from being accidentally withdrawn (symmetrical design allows tray insertion from either side)
- · Carerail fitted to either side support a wide range of accessories
- Translucent, impact-resistant trays & dividers, smooth surfaces & radiused edges assists cleaning (compatible with steam cleaning / high temperature machine washing)
- Trays feature an integral handle with a plastic clip-on label cover (coloured labels supplied white, blue, green, yellow & red)
- Each tray supplied with 2 x long & 2 x short dividers
- · Max load per tray 5kg
- Polyurethane ring buffers offer protection to both the trolley & fabric of the building
- 100mm pressed steel castors, 2 off braking

Dimensions (w x d x h)

- Overall Dimensions
 - Single Column 550 x 550 x 807mm
 - Double Column 1005 x 550 x 807mm
- Plastic Top
 - Single Column 455 x 455 (internal) at height 782mm
 - Double Column 915 x 455 (internal) at height 782mm
- Travs
 - Shallow Tray 355 x 405 x 100mm (internal)
 - Deep Tray 355 x 405 x 150mm (internal)

Standard Colours

Paint - Spearmint

Item	Description	Weight
CTM1/08/1S2D	1 x 100mm, 2 x 150mm Trays	20.9kg
CTM1/08/3D	3 x 150mm Trays	21.2kg
CTM1/08/3S1D	3 x 100mm, 1 x 150mm Trays	22.0kg
CTM1/08/4S	4 x 100mm Trays	21.6kg









CTM1/08/3S1D - shown with optional CT/WF, SBH/15 & CT/ALC

Options / Accessories

- Locking & Security page 54Overbridge & Tilt Bins page 60
- Defibrillator Platform page 61
- Carerail Mounted page 62

- Laptop Holders page 68
 Castors & Ring Buffers page 70
 Trays, Shelves, Dividers & Labels page 71



Caretray Trolleys - Mild Steel, Standard Level, Single & Double Column

Features

- Standard level, single & double column easy clean procedure trollevs
- · Quick & easy disassembly, assists thorough & deep cleaning
- Recyclable polypropylene one-piece top & side panels, smooth surfaces & radiused edges assists cleaning (compatible with steam cleaning)
 - Top Panel Raised edges on three sides provide defined grab points to act as push handles, whilst an integral rim at the front contains spillages (max load 30kg)
 - Side Panels Incorporate runners with positive stops preventing trays from being accidentally withdrawn (symmetrical design allows tray insertion from either side)
- Carerails fitted to either side support a wide range of accessories
- Translucent, impact-resistant trays & dividers, smooth surfaces & radiused edges assists cleaning (compatible with steam cleaning / high temperature machine washing)
- Trays feature an integral handle with a plastic clip-on label cover (coloured labels supplied white, blue, green, yellow & red)
- Each tray supplied with 2 x long & 2 x short dividers
- Max load per tray 5kg
- Polyurethane ring buffers offer protection to both the trolley & fabric of the building
- 100mm pressed steel castors. 2 off braking

Dimensions (w x d x h)

- Overall Dimensions
 - Single Column 550 x 550 x 993mm
 - Double Column 1005 x 550 x 993mm
- Plastic Top
 - Single Column 455 x 455 (internal) at height 968mm
 - Double Column 915 x 455 (internal) at height 968mm

- Shallow Tray 355 x 405 x 100mm (internal)
- Deep Tray 355 x 405 x 150mm (internal)

Standard Colours

• Paint - Spearmint







Options / Accessories Locking & Security - page 54 Overbridge & Tilt Bins - page 60 Defibrillator Platform - page 61 Carerail Mounted - page 62

- Laptop Holders page 68
 Castors & Ring Buffers page 70
 Trays, Shelves, Dividers & Labels page 71



Item	Description	Weight
CTM2/10/3S6D	3 x 100mm, 6 x 150mm Trays	46.0kg
CTM2/10/4S5D	4 x 100mm, 5 x 150mm Trays	45.6kg
CTM2/10/5S4D	5 x 100mm, 4 x 150mm Trays	45.3kg
CTM2/10/6S4D	6 x 100mm, 4 x 150mm Trays	46.7kg
CTM2/10/8D	8 x 150mm Trays	45.2kg
CTM2/10/8S2D	8 x 100mm, 2 x 150mm Trays	46.0kg



Caretray Trolleys - Mild Steel, High Level, Single & Double Column

Features

- High level, single & double column easy clean procedure trollevs
- · Quick & easy disassembly, assists thorough & deep cleaning
- Recyclable polypropylene one-piece top & side panels, smooth surfaces & radiused edges assists cleaning (compatible with steam cleaning)
 - Top Panel Raised edges on three sides provide defined grab points to act as push handles, whilst an integral rim at the front contains spillages (max load 30kg)
 - Side Panels Incorporate runners with positive stops preventing trays from being accidentally withdrawn (symmetrical design allows tray insertion from either side)
- · Carerail fitted to either side support a wide range of accessories
- Translucent, impact-resistant trays & dividers, smooth surfaces & radiused edges assists cleaning (compatible with steam cleaning / high temperature machine washing)
- Trays feature an integral handle with a plastic clip-on label cover (coloured labels supplied white, blue, green, yellow & red)
- Each tray supplied with 2 x long & 2 x short dividers
- Max load per tray 5kg
- Polyurethane ring buffers offer protection to both the trolley & fabric of the building
- 100mm pressed steel castors, 2 off braking

Dimensions (w x d x h)

- Overall Dimensions
 - Single Column 550 x 550 x 1410mm
 - Double Column 1005 x 550 x 1410mm
- Plastic Top
 - Single Column 455 x 455 (internal) at height 1385mm
 - Double Column 915 x 455 (internal) at height 1385mm
- Trays
 - Shallow Tray 355 x 405 x 100mm (internal)
 - Deep Tray 355 x 405 x 150mm (internal)

Standard Colours

· Paint - Spearmint

Item	Description	Weight
CTM1/14/2S5D	2 x 100mm, 5 x 150mm Trays	36.2kg
CTM1/14/3S4D	3 x 100mm, 4 x 150mm Trays	35.8kg
CTM1/14/4S3D	4 x 100mm, 3 x 150mm Trays	35.4kg
CTM1/14/5S3D	5 x 100mm, 3 x 150mm Trays	37.3kg
CTM1/14/6D	6 x 150mm Trays	35.0kg
CTM1/14/8S1D	8 x 100mm, 1 x 150mm Trays	37.2kg



- Options / Accessories

 Locking & Security page 54

 Carerail Mounted page 62

 Castors & Ring Buffers page 70

 Trays, Shelves, Dividers & Labels page 71



CTM2/14/10S6D - shown with optional CTD2/14SS (Double door fitted with Security Seal Lock)

Item	Description	Weight
CTM2/14/4S10D	2 x 100mm, 10 x 150mm Trays	60.6kg
CTM2/14/6S8D	6 x 100mm, 8 x 150mm Trays	62.8kg
CTM2/14/8S6D	8 x 100mm, 6 x 150mm Trays	62.3kg
CTM2/14/10S6D	10 x 100mm, 6 x 150mm Trays	65.0kg
CTM2/14/12D	12 x 150mm Trays	61.3kg
CTM2/14/16S2D	16 x 100mm, 2 x 150mm Trays	66.6kg



Caretray Racks - Mild Steel, Mobile, Single & Double Column

Features

- Mobile, single & double column easy clean racks
- Quick & easy disassembly, assists thorough & deep cleaning
- Recyclable polypropylene one-piece top & side panels, smooth surfaces & radiused edges assists cleaning (compatible with steam cleaning)
 - Top Panel Raised edges on three sides provide defined grab points to act as push handles, whilst an integral rim at the front contains spillages (max load 30kg)
 - Side Panels Incorporate runners with positive stops preventing trays from being accidentally withdrawn (symmetrical design allows tray insertion from either side)
- Translucent, impact-resistant trays & dividers, smooth surfaces & radiused edges assists cleaning (compatible with steam cleaning / high temperature machine washing)
- Trays feature an integral handle with a plastic clip-on label cover (coloured labels supplied white, blue, green, yellow & red)
- Each tray supplied with 2 x long & 2 x short dividers
- Max load per tray 5kg
- Polyurethane ring buffers offer protection to both the rack & fabric of the building
- 100mm pressed steel castors, 2 off braking

Dimensions (w x d x h)

- Overall Dimensions
 - Single Column 550 x 550 x 1782mm
 - Double Column 1005 x 550 x 1782mm
- Plastic Top
 - Single Column 455 x 455 (internal) at height 1757mm
 - Double Column 915 x 455 (internal) at height 1757mm
- Irays
 - Shallow Tray 355 x 405 x 100mm (internal)
 - Deep Tray 355 x 405 x 150mm (internal)

Standard Colours

• Paint - Spearmint

Item	Description	Weight
CRMM1/18/2S6D	2 x 100mm, 6 x 150mm Trays	41.2kg
CRMM1/18/3S6D	3 x 100mm, 6 x 150mm Trays	43.6kg
CRMM1/18/4S5D	4 x 100mm, 5 x 150mm Trays	42.3kg
CRMM1/18/5S5D	5 x 100mm, 5 x 150mm Trays	43.8kg
CRMM1/18/6S4D	6 x 100mm, 4 x 150mm Trays	43.4kg
CRMM1/18/8D	8 x 150mm Trays	42.8kg



- Options / Accessories

 Castors & Ring Buffers page 70

 Trays, Shelves, Dividers & Labels page 71



CRMM2/18/8S10D

Item	Description	Weight
CRMM2/18/4S12D	4 x 100mm, 12 x 150mm Trays	73.3kg
CRMM2/18/6S12D	6 x 100mm, 12 x 150mm Trays	76.3kg
CRMM2/18/8S10D	8 x 100mm, 10 x 150mm Trays	74.6kg
CRMM2/18/10S10D	10 x 100mm, 10 x 150mm Trays	77.6kg
CRMM2/18/12S8D	12 x 100mm, 8 x 150mm Trays	76.9kg
CRMM2/18/16D	16 x 150mm Trays	73.8kg



Caretray Racks - Mild Steel, Static, Single & Double Column

Features

- Mobile, single & double column easy clean racks
- Quick & easy disassembly, assists thorough & deep cleaning
- Recyclable polypropylene one-piece top & side panels, smooth surfaces & radiused edges assists cleaning (compatible with steam cleaning)
 - Side Panels Incorporate runners with positive stops preventing trays from being accidentally withdrawn (symmetrical design allows tray insertion from either side)
- Translucent, impact-resistant trays & dividers, smooth surfaces & radiused edges assists cleaning (compatible with steam cleaning / high temperature machine washing)
- Trays feature an integral handle with a plastic clip-on label cover (coloured labels supplied white, blue, green, yellow & red)
- Each tray supplied with 2 x long & 2 x short dividers
- Max load per tray 5kg
- Adjustable feet allow levelling on uneven floors

Dimensions (w x d x h)

- Overall Dimensions
 - Single Column 550 x 550 x 1670-1685mm
 - Double Column 1005 x 550 x 1670-1685mm
- Plastic Top
 - Single Column 455 x 455 (internal) at height 1645mm
 - Double Column 915 x 455 (internal) at height 1645mm
- - Shallow Tray 355 x 405 x 100mm (internal)
 - Deep Tray 355 x 405 x 150mm (internal)

Dimensions (w x d x h)

Paint - Spearmint

Item	Description	Weight
CRMS1/16/2S6D	2 x 100mm, 6 x 150mm Trays	38.4kg
CRMS1/16/3S6D	3 x 100mm, 6 x 150mm Trays	39.5kg
CRMS1/16/4S5D	4 x 100mm, 5 x 150mm Trays	39.5kg
CRMS1/16/5S5D	5 x 100mm, 5 x 150mm Trays	41.0kg
CRMS1/16/6S4D	6 x 100mm, 4 x 150mm Trays	40.6kg
CRMS1/16/8D	8 x 150mm Trays	39.6kg





Easy disassembly & cleaning

Options / Accessories
Trays, Shelves, Dividers & Labels - page 71





Item	Description	Weight
CRMS2/16/4S12D	4 x 100mm, 12 x 150mm Trays	70.2kg
CRMS2/16/6S12D	6 x 100mm, 12 x 150mm Trays	73.2kg
CRMS2/16/8S10D	8 x 100mm, 10 x 150mm Trays	72.3kg
CRMS2/16/10S10D	10 x 100mm, 10 x 150mm Trays	75.3kg
CRMS2/16/12S8D	12 x 100mm, 8 x 150mm Trays	74.6kg
CRMS2/16/16D	16 x 150mm Travs	71.5ka

Locking

- · Doors available as either
- assist cleaning)

- Four locking options
- Factory fitted



CAM Lock with Wing Handle

- Supplied with two keys
 Wing handle can operate independently of the key
 Locks keyed to differ, alike on request







CTM1/08/3S1D/CTD1/08CA

Description	Solid Grade Laminate	Clear Acrylic	Weight
Low Level, Single Column	CTD1/08CA	CTDC1/08CA	2.7kg
Low Level, Double Column	CTD2/08CA	-	5.3kg
Standard Level, Single Column	CTD1/10CA	CTDC1/10CA	3.5kg
Standard Level, Double Column	CTD2/10CA	-	6.8kg
High Level, Single Column	CTD1/14CA	CTDC1/14CA	5.3kg
High Level, Double Column	CTD2/14CA	-	10.3kg

CTM1/10/3S2D/CTDC1/10CA

Electronic Push Button Lock

- Uses a 10-button keypad
- Easy programming & code change
- Four operating codes with differing levels of authority Master / Sub-Master / User / Technician
- User codes are 4 digits
- Master & sub-master codes are 8 digits
- 15,000 openings on 2 x AAA cells (supplied)
- Low battery warning with failure override using 9v battery across LED's (+/-) & master code





CTM1/08/3S1D/CTD1/08EPB

Description	Solid Grade Laminate	Clear Acrylic	Weight
Low Level, Single Column	CTD1/08EPB	CTDC1/08EPB	2.9kg
Low Level, Double Column	CTD2/08EPB	-	5.5kg
Standard Level, Single Column	CTD1/10EPB	CTDC1/10EPB	3.7kg
Standard Level, Double Column	CTD2/10EPB	-	7.0kg
High Level, Single Column	CTD1/14EPB	CTDC1/14EPB	5.5kg
High Level, Double Column	CTD2/14EPB	-	10.5kg



CTM1/10/3S2D/CTDC1/10EPB

High Security Bolt-Lock

Features

Supplied with two keys







CTM1/10/3S2D/CTD1/10/HSB

Description	Solid Grade Laminate	Weight
Low Level, Single Column	CTD1/08HSB	2.9kg
Low Level, Double Column	CTD2/08HSB	5.6kg
Standard Level, Single Column	CTD1/10HSB	3.6kg
Standard Level, Double Column	CTD2/10HSB	7.1kg
High Level, Single Column	CTD1/14HSB	5.3kg
High Level, Double Column	CTD2/14HSB	10.6kg

Security Seal Lock

Features

• Accepts standard plastic security seals (order separately)







CTM1/08/3S1D/CTD1/08SS

Description	Solid Grade Laminate	Clear Acrylic	Weight
Low Level, Single Column	CTD1/08SS	CTDC1/08SS	2.6kg
Low Level, Double Column	CTD2/08SS		5.2kg
Standard Level, Single Column	CTD1/10SS	CTDC1/10SS	3.7kg
Standard Level, Double Column	CTD2/10SS		6.7kg
High Level, Single Column	CTD1/14SS	CTDC1/14SS	5.2kg
High Level, Double Column	CTD2/14SS		10.2kg

CTM1/10/3S2D/CTDC1/10SS

Security Clamp

Features

• Allows the trolley to be secured to the 'fabric of the building'



Important Note - Security Clamps

It is a requirement that when the trolley is not in use it should be secured to the 'fabric of the building'. Not only does this prevent casual theft, i.e. trolleys being pushed to an isolated area & broken into, but also encourages good nursing practice, giving the trolley a "permanent home".

"Safe & Secure Handling of Medicines - Royal Pharmaceutical Society, GB-March 2005"

Item	Description	Cable Length	Weight
MT/1/6F	Security Clamp	950mm	0.75kg

Security Seals

- Snap on plastic security seals
- Pack of 100



Item	Description
2PLCO026	Security Seals